Thermo Scientific Orion Star A211

Advanced pH benchtop meter

The Thermo Scientifio™ Orion™ Star A211 pH Benchtop Meter works for any lab or application. Its durability and dependability make it ideal for every sample, every day.









The Thermo Scientific Orion Star A211 pH Benchtop Meter is the ideal choice for any lab. Get the information you need quickly and easily from the large, backlit graphic LCD display. Informative display shows pH/mV/RmV and temperature with electrode status, time, date, sample ID, user ID and calibration points.

Features and Benefits

- Easy meter operation assured with on-screen text prompts, menu-specific function keys and a multi-language interface for English, Spanish, French, Italian, German and Chinese
- Don't miss a reading—AUTO-READ™
 locks in the stable reading on your screen,
 ready indicator alerts when readings are
 stable and timed reading gathers data in
 specific time intervals
- Up to 5 point pH calibration with automatic recognition for USA/NIST and DIN buffers or manual option for custom buffer values and calibration editing to fix individual points without a full recalibration

- Non-volatile memory holds up to 2000 data points with time and date stamp
- Easily transfer data and keep meter software up-to-date with the USB and RS232 ports and complimentary Orion Star Com computer software
- Stir samples with ease using the meter-controlled Orion Star stirrer probe (Cat. No. 096019)
- Included electrode arm and holder make it easier to maintain and place probes into samples
- Power the meter using the included universal power adapter or four AA batteries (sold separately)
- IP54-rated housing handles splashes and is wall-mountable to save space
- 3 year meter replacement warranty



Thermo Scientific Orion Star A211 pH Benchtop Meter

| | Specifications |
|-------------|--|
| рН | Range: -2.000 to 20.000 |
| | Resolution: 0.1, 0.01, 0.001 |
| | Relative Accuracy: ±0.002 |
| | Calibration Points: Up to 5 |
| | Calibration Editing: Yes |
| mV/RmV | Range: ±2000.0 mV |
| | Resolution: 0.1 mV |
| | Relative Accuracy: ± 0.2 mV or ± 0.05 % |
| | of reading whichever is greater |
| | E _H ORP Mode: Yes |
| Temperature | Range: -5 to 105 °C, 23 to 221 °F |
| | Resolution: 0.1 °C, 0.1 °F |
| | Relative Accuracy: ± 0.1 °C, ± 0.1 °F |
| | Offset Calibration: 1 point |
| Datalogging | Memory: 2000 with time and date stamp |
| | Log Function: Manual, ready |
| | (includes AUTO-READ), timed |
| Inputs | pH Electrode: BNC, reference pin |
| | ATC Probe: 8-pin mini-DIN |
| Output | RS232, USB |
| Power | AC Adapter: Included—universal, |
| | 100-240 VAC |
| | Battery Power: Optional—4 AAs |
| | Battery Life: 800 hrs |

| Meter Ordering Information | | |
|----------------------------|--|--|
| Part No. | Description | |
| STARA2110 | Orion Star A211 pH Benchtop Meter - Attached electrode stand and holder - Universal power adapter - Printed quick start guide and literature CD - USB computer cable | |
| STARA2115 | Orion Star A211 pH Benchtop Meter Kit - 8302BNUMD Orion ROSS Ultra Triode refillable, glass-body pH/ATC electrode - 810199 Orion ROSS pH buffer and solutions kit - Attached electrode stand and holder - Universal power adapter - Printed quick start guide and literature CD - USB computer cable | |

| STARA2116 | Orion Star A211 pH Benchtop Meter Kit - 8157BNUMD Orion ROSS Ultra Triode refillable, epoxy-body pH/ATC electrode - 810199 Orion ROSS pH buffer and solutions kit - Attached electrode stand and holder - Universal power adapter - Printed quick start guide and literature CD - USB computer cable |
|-----------|--|
| STARA2117 | Orion Star A211 pH Benchtop Meter Kit - 9157BNMD Orion Triode refillable, epoxy-body pH/ATC electrode - 910199 Orion pH buffer and solutions kit - Attached electrode stand and holder - Universal power adapter - Printed quick start guide and literature CD - USB computer cable |

| Ac | ccessories Ordering Information |
|----------------------|--|
| Part No. | Description |
| 096019 | Orion Star stirrer probe |
| STARA-BEA | Orion Star A series electrode arm, holder and benchtop meter-attachable base plate |
| STARA-HB | Orion Star A series freestanding base for electrode arm and holder |
| 8172BNWP | Orion ROSS Sure-Flow refillable, glass-body pH electrode |
| 8165BNWP | Orion ROSS Sure-Flow refillable, epoxy-body pH electrode |
| 8102BNUWP | Orion ROSS Ultra refillable, glass-body pH electrode |
| 927007MD | Orion stainless steel ATC temperature probe |
| 910104 | Orion pH 4.01 buffer, 475 mL |
| 910107 | Orion pH 7.00 buffer, 475 mL |
| 910110 | Orion pH 10.01 buffer, 475 mL |
| 810001 | Orion ROSS pH electrode storage solution, 475 mL |
| 910001 | Orion pH electrode storage solution, 475 mL |
| 9678BNWP | Orion Sure-Flow refillable epoxy-body ORP electrode |
| 967901 | Orion ORP standard solution, 475 mL |
| For more information | n, contact your local Thermo Scientific products dealer |

or call our customer and technical service experts at 1-800-225-1480 $\,$ (for the US and Canada) or visit www.thermoscientific.com/water.



Thermo Fisher Scientific Water Analysis Instruments Chelmsford, MA USA Quality Management System Registered to ISO 9001

thermoscientific.com/water

© 2013 Thermo Fisher Scientific Inc. All rights reserved. Certified System StandardsMark logo is a trademark of SAI Global. All other trademarks are the property of Thermo Fisher Scientific Inc. & its subsidiaries











Water Analysis Instruments

North America Toll Free: 1-800-225-1480 Tel: 1-978-232-6000 info.water@thermo.com

Netherlands

Tel: (31) 020-4936270 info.water.uk@thermo.com

Tel: (86) 21-68654588 wai.asia@thermofisher.com

India

Tel: (91) 22-4157-8800 wai.asia@thermofisher.com

Singapore Tel: (65) 6778-6876 wai.asia@thermofisher.com

Japan Tel: (81) 045-453-9175 wai.asia@thermofisher.com

Australia Tel: (613) 9757-4300 in Australia (1300) 735-295 InfoWaterAU@thermofisher.com



Redesigned electrode arm and holder

Part of Thermo Fisher Scientific